





On Dyna® Cervical Disk Protezi

*Metal on metal cervical disk protezi (Ti₆AlV₄)
Kanatsız cervical disk protezi*



- İç yüzeyler sürtünme katsayısını düşürecek TiN kaplama.
- Dış Yüzeyler Erken Kemik Tutunmasını sağlayacak Poroz Kaplama.
- Flexion, extantiona, lateral extantiona ve axial rotasyona müsaade eden dizayn
- Her yöne hareketliliğe müsaade eden dizayn
- Dislocke olmayı önleyen güvenli birleşim sistemi
- Gerektiğinde kolay çıkartmayı sağlayacak tutucu kanallar.
- Tek seferde implantasyon.

66.550.006	On DYNA Cervical Disk Proshesis	Small (5 mm)
66.550.008	On DYNA Cervical Disk Proshesis	Medium (6 mm)
66.550.010	On DYNA Cervical Disk Proshesis	Large (7 mm)
66.550.012	On DYNA Cervical Disk Proshesis	X Large (8 mm)
66.550.014	On DYNA Cervical Disk Proshesis	XX Large (9 mm)



TiN Kaplama Nedir?

Coating Name	Coating	Hardness (HV 0.05) a	Adhesion (N) b	Oxidation Temperature (°C/°F)	Friction Coefficient c	Thickness	Surface Roughness (R μM) d	Color
Medthin-Gold	TiN	2900±200	70	500°C 950°F	0.6	2 to 4	0.2	Gold 