



CE
0123

REF A.5001

QUADRIJET®

TEKNİK TABLO

Perianal anestezi için tek kullanımlık dörtlü luer iğne - şırınga bağlantısı

Quadrijet tek kullanımlık iğne şırınga bağlantısı "minör" dar alanlı proktolojik cerrahi müdahalelerde perianal anesteziler için geliştirilmiştir. Bu teknik dört ayrı yerden aynı şırıngayla yapılan anestezinin yerini almıştır.

Dört iç kanalın ölçüleri dizaynı anestezik çözeltinin dört ayrı bölgeye aynı alanda enjekte edilmesini sağlayarak, anestezinin hasta için daha dayanılabilir olmasını sağlar.

ETO ile steril edilmiştir.

Kutu Adedi: 100 Adet

TECNICAL SHEET

**Quadruple luer needle-syringe
connector for perianal anaesthesia.
Sterile. Disposable.**

The Quadrijet disposable quadruple luer needle-syringe connector has been developed for perianal anaesthesia in minor proctologic operations. This technique replaces the need to carry out separate injections at four sites, with one single injection.

The gauged design of the four internal channels means that the anaesthetic solution is injected in even doses in each quadrant, so making anaesthesia more tolerable for the patient.

Sterilised by ETO

Quantity per box: 100x pieces

QUADRIJET®

CE
0123



Perianal anestezi için
dörtlü luer iğne - şırınga
bağlantısı.
Steril & Tek kullanımlık.

- Quadruple luer
- needle-syringe for
- connector perianal
- anaesthesia.
- Sterile & Disposable.

